


## **ECG ELECTRODES NEONATE / PEDIATRIC**



- ***Suitable to the current needs and specific conditions***

As the bio-adhesive gel ensures the adhesive properties needed in the incubator humidity, the electrodes do not require a frequent replacement. The electrodes and cable are very soft and do not cause any irritation. The peaceful rest of the baby is not disturbed, your work neither.

- ***Suitable for MRI and X-Ray procedures***

Intégral Process electrodes are radio translucent and do not show any shadow on the MRI screen or X-Ray.

- ***Easy to use***

The electrodes are in pack of 3 (exact number for one patient). They are easy to unpack and place. The gel does not stick to clothes or sheets.

- ***Display box***

The display box allows an optimal storage and avoids any packs mixing. The packs are also very easy to take from the box.



## FEATURES :


- Latex free electrodes
- Low impedance
- One size for all body types
- PVC free electrodes
- Phtalate free wire

## MATERIAL :

- Biocompatibility
- Contact Ag/AgCl (Silver/Silver chloride)
- Solid gel can be repositioned
- Radio translucent wire
- Polyester cloth 25 X 22 mm
- Contact type : nickel plated brass
- Conductor wire : carbon
- Coloured polyethylene coating

## PACKAGING :

- Pouch of 3 electrodes stripe
- 50 pouches per box  
(box of 150 electrodes)




## Part N°40556IP



**4MM SOCKET**

**Length 0.60m - Radiotranslucent**

## Part N°40554IP



**DIN 1,5 mm plug**

**Length 0.60m - Radiotranslucent**

## Part N°40557IP



**DIN 1,5 mm plug**

**Length 0.50m - Copper**