

FOETAL TRANSDUCERS

Make the difference !



A single source to answer all your needs in terms of sensors

- Our sensors are CE (CE0086) and FDA approved. They allow ensuring a reliable compatibility.
- Our sensors are waterproof.
- The flat shape has been designed to avoid any hurt to the patient.
- Intégral Process offers its expertise for the repair of your defective sensors out of warranty.
- Our sensors have a 1 year warranty except for Philips Avalon.
- A wide range covering most of the foetal monitors.





















- **PHILIPS old range**
- **PHILIPS Avalon**
- **G.E. COROMETRICS**
- **SONICAÏD**



www.integral-process.com

Capter la vie - Sensing Life

PHILIPS (old generation) – Compatible sensors			
STRAP SIDE	PART NUMBER / TYPE	CONNECTOR	PATIENT SIDE
	<p>TOCO FM10841IP Original part N° M1355A Serie 1350 – Viridia 50 FX</p>	 Brown connector	
	<p>ULTRASOUND FM10835IP Original part N° 1356A Serie 1350 – Viridia 50 FX</p>	 Pink connector	
	<p>TOCO ATG100-T Original part N° 15248 Serie 8040</p>	 Brown connector	
	<p>ULTRASOUND ATG100-U Original part N° 15245 Serie 8040</p>	 Blue connector	
PHILIPS (new generation) – Original sensors			
STRAP SIDE	PART NUMBER / TYPE	CONNECTOR	PATIENT SIDE
	<p>ULTRASOUND E164CUA Original part N° M2736A AVALON FM 20-30-40-50</p>	 Pink connector	
	<p>TOCO E164CTA Original part N° M2734A AVALON FM 20-30-40-50</p>	 Pink connector	
	<p>TOCO + E164CTB Original part N° M2735A AVALON FM 20-30-40-50</p>	 Pink connector	

SONICAID - Compatible sensors			
STRAP SIDE	PART NUMBER / TYPE	CONNECTOR	PATIENT SIDE
	<p>ULTRASOUND E235AUB Original part N° 8400-6919 OXFORD 1,5 MHz</p>	 Yellow connector	
	<p>ULTRASOUND E235AUA Original part N° 8400-6920 OXFORD 2 MHz</p>	 Blue connector	
	<p>TOCO E235ATA Original part N° 8400-6920 TEAMDUO</p>	 Pink connector	
GE / COROMETRICS - Compatible sensors			
STRAP SIDE	PART NUMBER / TYPE	CONNECTOR	PATIENT SIDE
	<p>ULTRASOUND FM10833IP Original part N° 5700 LAX NAUTILUS 116/118/120/170/259/340T</p>	 Blue connector	
	<p>ULTRASOUND FM10834IP Original part N° 5700 HAX NAUTILUS 116/118/120/170/259/340T</p>	 Blue connector	
	<p>TOCO FM10839IP Original part N° 2264 LAX NAUTILUS 115/116/118/120/145/150 151/155/170/259/340T</p>	 White connector	
	<p>TOCO FM10840IP Original part N° 2264 HAX NAUTILUS 115/116/118/120/145/150 151/155/170/259/340T</p>	 White connector	