



**NEW SINGLE PATIENT NEONATE SPO<sub>2</sub>  
SENSOR**

**Can be repositioned thanks to its Velcro fastening**



[www.integral-process.com](http://www.integral-process.com)

Capter la vie - Sensing Life

# DisposableWrap – The new SPO<sup>2</sup> disposable sensor

EnviteC, Leader in the technology of non-invasive medical sensors offers a new disposable SPO<sup>2</sup> sensor **DisposableWrap**.

## Main Features :

- Patient safety
- Allergy proof
- Improvement for placement and repositioning on the patient
- Increases the patient comfort



More than all other EnviteC sensors, **DisposableWrap** SPO<sup>2</sup>sensor allows a great comfort during the repositioning with patient safety.

## The best safety for the patient thanks to a high quality comfort :

- Hypoallergenic Velcro® wrap agreed for medical application allowing an optimal comfort.
- Latex free
- Skin friendly thanks to a very soft material
- Reduces the sweating
- High quality components
- Avoids injuries during repositioning.
- Ergonomic design for an easy placement.

## Optimal handling and placement easiness for the medical staff.

- Easy to place and reposition
- No storage availability limit

## Technical Datas

Sensor	Neonate
Cable	90 cm

Compatibility Chart	Sensors part numbers	Adapter cables
BCI	92880ME	90742BCI
GE Ohmeda	92280ME	90188IP
GE connecteur rond	92280ME	90190IP
GE connecteur rectangulaire	92280ME	90187IP
Nellcor® *	92780ME	*
Nihon Kohden	92680ME	90775
Novamatrix	92580ME	90748IP
Ohmeda	92480ME	90738IP

\* Only for Nellcor® type technologies

Manufactured by Envitec Honeywell